

WELLS
JOHNSON
COMPANY

Stat Kit® 900

Part# 20-5080

The Stat Kit 900 (Item #900) is priced at \$1250.00, and is the most extensively equipped Stat Kit. In addition to the multitude of drugs, the Stat Kit 900 contains emergency oxygen, and Ambu disposable resuscitator, and a manual suction unit. You have the drugs, oxygen, and equipment to start treatment on virtually any life-threatening emergency associated with an office-based practice.

The Stat Kit 900 is portable at 28 pounds. With our newly redesigned foam tray protection, all drug contents are cushioned snugly in place.

You might also want to consider the optional ACLS drug package (Item# 20-5082). This package contains two vials each of Amiodarone and Vasopressin.

The oxygen system has a preset flow rate for speed of setup and administration. It is preset to deliver a minimum of 8LPM for approximately 30 minutes. The medical "C" spun aluminum cylinder is easily refilled at oxygen supply houses.

The Ambu SPUR disposable resuscitator is one of the most popular selling disposable manual resuscitators. There is no sterilization and no difficulty in assembly or loss of components. It connects easily to the oxygen system.

The manual suction unit included in the Stat Kit 900 is compact, with a reusable pump handle. The disposable container and catheter traps more than 200ml, quickly suctioning fluids from the nose and mouth. The disposable container is easily sealed and discarded.



8000 S Kolb Rd, Tucson AZ 85756

phone 800.528.1597

fax 520.885.1189

sales@wellsgrp.com

www.wellsgrp.com

Stat Kit® 900

Part# 20-5080

Dimensions: 6 3/8" x 23" x 27"

Weight: 28 lbs

CONTENTS OF STAT KIT 900

Drugs-Prefilled Syringes:

Atropine (2)
Dextrose (1)
*Diazepam (2)
Epinephrine 1 1/2" needle (3)
Lidocaine (2)
Sodium Bicarbonate (1)

Drugs-Ampules & Vials:

Diphenhydramine (2)
Epinephrine (3)
Furosemide (2)
Lanoxin (1)
Nalbuphine (2)
Naloxone (2)
Procainamide (1)
Promethazine (2)
Solu-Cortef (1)

Drugs-Miscellaneous:

Activated Charcoal (1)
Albuterol Inhaler (1)
Ammonia Inhalants (2)
Aspirin Tablets (4)
Clonidine Tablets (2)
Diphenhydramine Caplets (4)
Nitroglycerin Tablets (25)

IV Equipment:

0.9% Sodium Chloride 500ml (1)
I.V. Administration Tubing (1)
I.V. Catheter 20ga x 1" (2)

Safety Syringes & Needles:

10cc w/21ga x 1 1/2" (1)
10cc w/22ga x 1 1/2" (1)
3cc w/20ga x 1 1/2" (1)
3cc w/21ga x 1 1/2" (1)
3cc w/22ga x 1" (1)
3cc w/25ga x 5/8" (1)

Resuscitation Equipment:

Disposable Resuscitator
w/ Adult Mask (1)
Endotracheal Tubes
3mm w/stylet (1)
5mm w/stylet (1)
7mm w/stylet (1)
Laryngoscope, Large (1)
Laryngoscope, Small (1)
(Interchangeable light source)
Manual Suction Unit (1)
Oropharyngeal Airways
Small (1)
Medium (1)
Large (1)

Oxygen Equipment:

Oxygen Mask, with strap and
non-rebreathing bag (1)
Oxygen Cylinder with
toggle handle (1)
Roll of Tubing (1)

Surgical Instruments:

Alcohol Sponges (2)
Gauze Sponges 3x3 (2)
Gloves, Nitrile (2)
Hemostat (2)
Needle Holder (1)
Scalpel (1)
Scissors (1)
Sutures
4-0 Prolene (1)
2-0 Vicryl (1)
Tape (1)
Thumb Forceps (1)
Tourniquet (1)

Monitoring Equipment:

Sphygmomanometer (1)
Stethoscope (1)

Organizational Materials:

AHA Algorithm Book (1)
Plastic Seals (4)
Quarterly Inspection Record (1)
Reorder Form (1)
Treatment Tag (1)

Optional ACLS Package (20-5082)

Amiodarone (2)
Vasopressin (2)

**DEA certificate required for shipment of controlled substance.*