

ELITE Infiltrator/Aspirator

SYSTEM INCLUDES:

- Power Cord
- Canister Holder Ring
- Canister Adapter
- Aspirator Foot Switch
- Pneumatic Infusion Foot Switch
- IV Pole
- 10 pcs Sterile, Disposable Suction Tubing
- 10 pcs Sterile, Disposable Infusion Tubing
- 10 pcs Disposable Canister Liners
- 2 Bacterial Filters (1 installed in machine)
- 2 Infiltration Handles
- 1 Overflow Trap Assembly (Pre-Installed)
- 2 Infiltration Cannulas



Part# 20-5061-00 110V
20-5066-00 220V

The Elite combines high powered infiltration and aspiration into a single vertical unit with a small footprint.

- High Infusion Rate 0-1200 ml/min
- Precise Settings and Controls
- High Flow Rate 6.8 cubic ft/min
- Eye Level Controls
- Small Footprint 18"x18"

Infiltrator Specifications

- Flow Rate: 0-1200 ml/min
- Pressure: 25 PSI continuous
40 PSI intermittent

Aspirator Specifications

- Maximum Vacuum: 29+ in Hg
- Maximum Flow Rate: 6.8 cubic ft/min
- Number of Pumps: 3
- Horsepower: 1HP

ACCESSORY KIT OPTIONS

BAXTER™

1500CC SYSTEM

Baxter™ Canister
Reusable 1500cc
20-5153-00

Baxter™ Liner
Disposable 1500cc
20-5154-00 (Case of 10)

ASPIRATION TUBING

Standard Ultra-Limp™
Length 9ft
24-5103-00 (Case of 10)

Flex-Tube™ B
Lightweight, 12 ft
24-5102-02 (Case of 10)

CANISTER HOLDERS

Baxter Canister Holder fits
on top panel of the aspirator
to hold 1500cc Baxter Canisters in place.
20-5159-01

SORENSON™

1900CC SYSTEM

Sorenson™ Canister
Reusable 1900cc
20-5157-00

Sorenson™ Liner
Disposable 1900cc
20-5158-00 (Case of 10)

ASPIRATION TUBING

MB Ultra-Limp™
Length 9ft
24-5104-00 (Case of 10)

CANISTER HOLDERS

Sorenson Canister Holders fits
on top panel of the aspirator
to hold 1900cc Sorenson™
Canisters in place.
20-5159-00

QUICK-FIT™

3000CC SYSTEM

Quick-Fit™ Canister
Reusable 3000cc
20-5156-00

Quick-Fit™ Liner
Disposable 3000cc
20-5155-00 (Case of 10)

ASPIRATION TUBING

MB Ultra-Limp™
Length 9ft
24-5104-00 (Case of 10)

CANISTER HOLDERS

Quick-Fit Canister Holders fits
on top panel of the aspirator
to hold 3000cc Quick-Fit™
Canisters in place.
20-5159-03

HERCULES

Our top-of-the-line aspirator has unprecedented performance and quality. HERCULES™ features state-of-the-art technology in exterior component design with polymer casing that's scratch-resistant and defies denting. HERCULES™, with three high-powered dual-cylinder piston pumps, is the most powerful aspirator on the market and yet its operation is quiet. Standard all brass fittings resist corrosion and pre-tested high pressure tubing exceeds all industry requirements ensuring reliability for years to come.

- Most powerful aspirator available
- Dent-free scratch resistant casing
- Superior quality inside and out
- Available in 110V or 220V
- 5-year Warranty on parts and labor



PART #: 20-5007-00 (110V)
PART #: 20-5010-00 (220V)

ITEM NUMBER: 20-5007-00 (110V)
20-5010-00 (220V)

- Weight 73 lbs
- Dimensions 18 1/2" 18" x 15"
- Maximum Vacuum 29+ inches Hg
- Maximum Flow Rate 6.8 cubic feet per minute
- Sound Level 53+/-3 db @ 6ft
- Pump type 2 cylinder piston
- Number of pumps 3
- Horsepower 1 hp

PARTS & ACCESSORIES

- Foot Switch (single) 20-5102-00
- Dual Footswitch 20-5103-00
- Bacterial Filter 20-5201-00
- Overflow Trap (plastic) with:
6" and 9" silicone tubing 17-5100-00
- 6" Silicone Tubing 17-5101-00
- 9" Silicone Tubing 17-5106-00

ACCESSORY KIT OPTIONS

BAXTER™

1500CC SYSTEM

Baxter™ Canister
Reusable 1500cc
20-5153-00

Baxter™ Liner
Disposable 1500cc
20-5154-00 (Case of 10)

ASPIRATION TUBING

Standard Ultra-Limp™
Length 9ft
24-5103-00 (Case of 10)

Flex-Tube™ B
Lightweight, 12 ft
24-5102-02 (Case of 10)

CANISTER HOLDERS

Baxter Canister Holder fits
on top panel of the aspirator
to hold 1500cc Baxter Canisters in place.
20-5159-01

SORENSON™

1900CC SYSTEM

Sorenson™ Canister
Reusable 1900cc
20-5157-00

Sorenson™ Liner
Disposable 1900cc
20-5158-00 (Case of 10)

ASPIRATION TUBING

MB Ultra-Limp™
Length 9ft
24-5104-00 (Case of 10)

CANISTER HOLDERS

Sorenson Canister Holders fits
on top panel of the aspirator
to hold 1900cc Sorenson™
Canisters in place.
20-5159-00

QUICK-FIT™

3000CC SYSTEM

Quick-Fit™ Canister
Reusable 3000cc
20-5156-00

Quick-Fit™ Liner
Disposable 3000cc
20-5155-00 (Case of 10)

ASPIRATION TUBING

MB Ultra-Limp™
Length 9ft
24-5104-00 (Case of 10)

CANISTER HOLDERS

Quick-Fit Canister Holders fits
on top panel of the aspirator
to hold 3000cc Quick-Fit™
Canisters in place.
20-5159-03

WHISPERATOR

The WHISPERATOR™, featuring a smooth grey casing, manufactured in dent-free, scratch resistant polymer plastic. The interior has been re-designed with simplicity in mind: powered by dual-synchronized piston pumps and surrounded by a special 2" sound absorbing foam. The WHISPERATOR™ delivers time-saving performance without being heard. It's truly the quietest aspirator Wells Johnson has ever manufactured. When you choose the WHISPERATOR™, you'll receive superior quality throughout: all brass fittings and pre-tested high-pressure tubing, exceeding all industry requirements and ensuring reliability for years to come.



PART#: 20-5026-00 (110V)
PART#: 20-5028-00 (220V)

- ITEM NUMBER:** 20-5026-00 (110V)
 20-5028-00 (220V)
- Weight 64 lbs
 - Dimensions 18 1/2" x 18" x 15"
 - Maximum Vacuum 29+ inches Hg
 - Maximum Flow Rate 5.2 cubic feet per minute
 - Sound Level 50+/-3db @ 6ft
 - Pump Type 2 cylinder piston
 - Number of pumps 2
 - Horsepower 2/3 hp

PARTS & ACCESSORIES

- Foot Switch (single) 20-5102-00
- Dual Footswitch 20-5103-00
- Bacterial Filter 20-5201-00
- Overflow Trap (plastic) with:
 6" and 9" silicone tubing 17-5100-00
- 6" Silicone Tubing 17-5101-00
- 9" Silicone Tubing 17-5106-00

ACCESSORY KIT OPTIONS

BAXTER™

1500CC SYSTEM

Baxter™ Canister
 Reusable 1500cc
 20-5153-00

Baxter™ Liner
 Disposable 1500cc
 20-5154-00 (Case of 10)

ASPIRATION TUBING

Standard Ultra-Limp™
 Length 9ft
 24-5103-00 (Case of 10)

Flex-Tube™ B
 Lightweight, 12 ft
 24-5102-02 (Case of 10)

CANISTER HOLDERS

Baxter Canister Holder fits on top panel of the aspirator to hold 1500cc Baxter Canisters in place.
 20-5159-01

SORENSEN™

1900CC SYSTEM

Sorenson™ Canister
 Reusable 1900cc
 20-5157-00

Sorenson™ Liner
 Disposable 1900cc
 20-5158-00 (Case of 10)

ASPIRATION TUBING

MB Ultra-Limp™
 Length 9ft
 24-5104-00 (Case of 10)

CANISTER HOLDERS

Sorenson Canister Holders fits on top panel of the aspirator to hold 1900cc Sorenson™ Canisters in place.
 20-5159-00

QUICK-FIT™

3000CC SYSTEM

Quick-Fit™ Canister
 Reusable 3000cc
 20-5156-00

Quick-Fit™ Liner
 Disposable 3000cc
 20-5155-00 (Case of 10)

ASPIRATION TUBING

MB Ultra-Limp™
 Length 9ft
 24-5104-00 (Case of 10)

CANISTER HOLDERS

Quick-Fit Canister Holders fits on top panel of the aspirator to hold 3000cc Quick-Fit™ Canisters in place.
 20-5159-03

GENERAL ASPIRATOR

Suited for smaller procedures.
Very portable and lightweight.
Maintenance Free.

ITEM NUMBER: 20-5000-00 (110V)
20-5001-00 (220V)

- Weight 36 lbs
- Dimensions 14.5" x 9" x 16"
- Maximum Vacuum 29+ inches Hg
- Maximum Flow Rate 2.0 cubic feet per minute
- Sound Level 48+/-3db
- Pump Type Two-cylinder piston
- Number of pumps 1
- Horsepower 1/3 hp



PART #: 20-5000-00 (110V)
PART #: 20-5001-00 (220V)

ACCESSORY KIT OPTIONS

BAXTER™

1500CC SYSTEM

Baxter™ Canister
Reusable 1500cc
20-5153-00

Baxter™ Liner
Disposable 1500cc
20-5154-00 (Case of 10)

ASPIRATION TUBING

Standard Ultra-Limp™
Length 9ft
24-5103-00 (Case of 10)

Flex-Tube™ B
Lightweight, 12 ft
24-5102-02 (Case of 10)

CANISTER HOLDERS

Baxter Canister Holder fits on top panel of the aspirator to hold 1500cc Baxter Canisters in place.
20-5159-01

SORENSON™

1900CC SYSTEM

Sorenson™ Canister
Reusable 1900cc
20-5157-00

Sorenson™ Liner
Disposable 1900cc
20-5158-00 (Case of 10)

ASPIRATION TUBING

MB Ultra-Limp™
Length 9ft
24-5104-00 (Case of 10)

Flex-Tube™ A
Lightweight 12 ft
24-5102-01 (Case of 10)

CANISTER HOLDERS

Sorenson Canister Holders fits on top panel of the aspirator to hold 1900cc Sorenson™ Canisters in place.
20-5159-00

QUICK-FIT™

3000CC SYSTEM

Quick-Fit™ Canister
Reusable 3000cc
20-5156-00

Quick-Fit™ Liner
Disposable 3000cc
20-5155-00 (Case of 10)

ASPIRATION TUBING

MB Ultra-Limp™
Length 9ft
24-5104-00 (Case of 10)

Flex-Tube A™
Lightweight, 12 ft
24-5102-01 (Case of 10)

CANISTER HOLDERS

Quick-Fit Canister Holders fits on top panel of the aspirator to hold 3000cc Quick-Fit™ Canisters in place.
20-5159-03